
8D HIFU

High Intensity
Focused Ultrasound

HIFU

HIFU, a non-invasive treatment, goes deeper to stimulate production of your body's own collagen and elastin, for natural-looking results without surgery or downtime.


3 = 10 years younger
Session




Theory

What is 8D HIFU?

High Intensity Focused Ultrasound (HIFU) delivers the heat energy to SMAS layers and subcutaneous tissue that can stimulate and renew the skin's collagen and therefore improving the surface and reducing sagging of the skin. This technique can be applied to the face as well as the whole body, and also, it works equally well on every skin tone.




01



Before the treatment, collagen in our skin is weak and disorganized which makes our skin saggy with fine lines or wrinkles.

02



Micro focused ultrasonic rays or energy passes through the skin tissue without damaging the cells in skin layer. It stimulates the collagen in skin.

03






The skin responds to the stimulations, make up new collagen and elastin at the same time, so the skin becomes more tightening and recover youthful state after one month.

Handpiece & Interface

01 5 Cartridge

5 cartridge can penetrate different depth of face and body precisely, multifunctional, a real low invest, high return machine.



Handpiece & Interface

02 Friendly Software



2 handles systems, clear to choose



professional system, easy to set parameter



8D system, shown exact treatment area

Application

Suitable for what kind of conditions ?

Treat sagging skin, aging phenomenon for any areas such as forehead, eye brow, cheekbones, apples muscle, cheeks, chin, neck and body



Advantage

Quick & short

Only 30mins one face treatment

SMAS contraction

Collagen remodeling, elastine fiber contraction

Lasting effect

Good result will last 2-3 years at least

None downtime

Skin just become red within the first several hours, then skin recover

Instant result

After treatment you will see 20-30% immediate result, then 3-6 months achieve best result

Totally non-invasive



Benefits



HIFU can generate up to
\$500/session

Treatment site	Average Fee	Procedures	New Revenue
Face	\$300	2	\$600
Arm	\$200	2	\$400
Abdomen	\$300	2	\$600
Thigh	\$300	2	\$600
Calf	\$200	2	\$400

Testimonials



Dr Justin Boey Medical Director

Beauty Clinic Owner

My clients like the treatment very much. They enjoy the process very much. I am very happy to see my clients skin get better.



Dr Haus Ranu

Beauty Clinic Owner

Our clinic use 8D HIFU machine for 3 years, it have the same result with Ultherapy, obvious lifting result only after one time treatment, HIFU help us attracted more and more cut-somers, love it!



Norman Bethune

Beauty Salon Owner

The results are simply incredible specially in the area around my cheeks and mouth where I can see a constant smoothness and disappearance of fine and expression lines. Also the wrinkles around my eyes and forehead have been substantially minimised. You see results in 3 weeks with Botox.

Before & After



Technical Specification

Irradiation Type	High Intensity Focused Ultrasound
Frequency	4MHz, 7MHz,
Max Energy	2.0J/m ²
LCD Screen	15 Inch True Color Touch LCD ADSS System
Standard Cartridges for single line handle	4MHz: 4.5mm; (cheek/neck)
	7MHz 3mm(forehead/cheek/neck)
Standard Cartridges for 12 lines handle	3MHz 7mm(body) 2MHz 13mm(body)
Optional cartridges for single line handle	2MHz 13mm; (body) 3MHz 7mm;(body)
	10MHz 1.5mm(eyes area)
Optional cartridges for 12 lines handle	4MHz 4.5mm(body) 7MHz 3mm(body)
Demensions of Machine	45*30*15cm
N.W./G.W.	20/22KG
Voltage	AC110-220V 50/60Hz