

WHY ACTION II?

- Variety of aesthetic indications and gynecological approaches
- Total Fractional Capabilities
- Long pulse mode: Pulse duration of 1,000 ms
- Shining Peel for Gentle Rejuvenation

ACTION II GUI



SYSTEM SPECIFICATIONS

Laser Type	Er:YAG 2940 nm	
Operating Mode	Dermatology	Normal Mode LongPulse Mode Fractional Mode
	Gynecology (Petit Lady)	Multiple Mode LongPulse Mode Dual Mode
Beam Profile	Top Hat	
Pulse Width	Normal Mode: 250 μs LongPulse Mode: max 1000 ms	
	Fractional H/P	12 mmØ, 9x9 mm
Handpiece	Zoom H/P (Optional)	1 – 7 mmØ
	Shining Peel H/P (Optional)	14 mmØ
	Petit H/P	360° 90°
Handpiece Auto Sensing	Yes	
Pulse Rate	Single, 1, 2, 5 and 10 Hz	
User Interface	8 inch LCD Display / Touch screen	
Memory	*Shining Peel* & 8 point	
Aiming Beam	Diode 655 nm (Red), 1 mW, Adjustable Brightness	
Dimensions (mm) / Without Arm	872 (H) x 656 (L) x 295 (W)	



ACTION II
 4 in 1 Innovation
 of Er:YAG Aesthetic Laser
 DUAL FRACTIONAL
 REJUVENATION & SHINING

LUTRONIC WORLD CENTER
 219 Sowon-ro, Deogyang-gu, Goyang-si, Gyeonggi-do, Korea
 TEL: +82-31-908-3440 / FAX: +82-1600-3660
 E-MAIL: korea@lutronic.com

LUTRONIC GLOBAL CENTER
 19 Fortune Drive, Billerica, MA 01821 USA
 TEL: +1-888-588-7644
 E-MAIL: office@lutronic.com / usa@lutronic.com

LUTRONIC GERMANY
 Esplanade 41, 20354 Hamburg, Germany
 TEL: +49 (0)40 6963 999 - 50 / FAX: +49 (0)40 6963 999 - 59
 E-MAIL: germany@lutronic.com

LUTRONIC FRANCE
 23 Rue Francois Jacob, 92500 Rueil-Malmaison, France
 TEL: +331 47 51 12 10 / E-MAIL: france@lutronic.com

LUTRONIC JAPAN
 12F, Tokyo Front Terrace, 2-3-14, Higashi-Shinagawa, Shinagawa-ku, Tokyo, Japan
 TEL: +81-3-6433-3041 / FAX: +81-3-6433-3042 / E-MAIL: japan@lutronic.com

LUTRONIC CHINA
 Yuzhou plaza 2001 room, Jingang road 299 no, Pudong district, Shanghai, 200120, China
 TEL: +86-21-5019-3891,3893 / FAX: +86-21-5019-3865
 E-MAIL: china@lutronic.com

©2019, LUTRONIC Group of Companies. All rights reserved. LUTRONIC, its logo, ACTION II are trademarks or registered trademarks of the LUTRONIC Group of Companies. *Specifications are subject to change without notice. This product or the use of this product is covered by one or more US and foreign patents or patent application pending. ACTION II BR 20190830
 For International Use Only

