

power to practice

HS-900 Multi-Technology Laser Platform

THE COMPLETE 8-IN-1 LASER SYSTEM FOR MODERN CLINICS

Exclusively Distributed in Australia by Power to Practice





Built to grow your practice

The HS-900 is a powerful and adaptable laser platform that combines eight energy-based technologies in one system. Designed to replace the need for multiple standalone machines, it enables aesthetic clinics to deliver a wide range of treatments using a single device — with seamless transitions between handpieces and a highly intuitive user interface. Built for clinics that demand versatility, clinical precision, and long-term value.

Key features include:

- 8 interchangeable laser and energy-based handpieces in one system
- Automatic handpiece detection for seamless operation
- Powerful energy output across all modalities
- TGA-listed, CE-certified, and FDA-cleared
- Smart, upgradeable software platform
- High satisfaction rates from patients and clinical users

Specifications

Model NO.: HS-900
Brand Name: APOLOMED
OEM/ODM: Professional Design Team and Rich Manufacture Experience
Certificate: ISO 13485, SGS ROHS, TUV Medical CE, TGA



Interface: 8" true colour touch screen \$24,990 AUD
Dimensions: 65*48*115cm (L*W*H)
Weight: 72kg
\$24,990



2940NM ER:YAG FRACTIONAL ABLATIVE LASER \$2,890 AUD
Spot size: 10x10mm, Φ 6mm, Φ 9mm, Φ 1-3.4mm
Pulse width: 1.5ms, 2ms, 2.5ms
Energy: 7x7 pixel: 12-28mJ/MTZ
9x9 pixel: 4-20mJ/MTZ
Beam expander: 500-2000mJ
Zoom lens: 500-2000mJ
Application: Skin resurfacing, wrinkles and fine lines, photodamage, texture irregularities, wart and nevus removal.



1064NM LONG PULSE ND:YAG LASER \$1,990 AUD
Pulse width: 10-40ms
Repetition rate: 0.5-1Hz
Energy: Φ 9mm: 10-110J/cm² Φ 6mm: 60-260J/cm² Φ 2.2*5mm: 150-500J/cm²
Application: Hair removal for all skin types, leg veins, vascular lesions, and wrinkle removal.



1064/532NM Q-SWITCH ND:YAG LASER

\$1,990 AUD

Spot size: 1-5mm
Pulse width: <10ns (single pulse)
Repetition rate: 1-10Hz
Maximum energy: 2400mJ-4700mJ
Application: Tattoo and tattoo wound removal, eyebrow and lip line removal, pigmented epidermal and dermal lesions, vascular lesions (telangiectasia), and soft peeling.



1540NM ER:GLASS FRACTIONAL LASER

\$1,990 AUD

Pulse width: 0ms, 15ms
Energy: 10×10 pixel: 30-70mJ/MTZ
18×18 pixel: 5-25mJ/MTZ
Application: Skin resurfacing, surgical scars, acne scars, stretch marks, melasma, and wrinkles.



IPL SHR/EPL

\$1,490 AUD

Spot size: 15*50mm
Wave length: 420-1200nm
Filter: 420/510/560/610/640-1200nm, SHR filter
Energy: 10-60 level
Application: Permanent hair removal, acne removal, skin toning, skin rejuvenation, epidermal pigment removal, spot and freckle removal, skin firming and tightening, and vascular treatment.



RF MONOPOLAR OR RF BIPOLAR

\$1,490 AUD

Output power: 200W
RF tip: Φ18mm, Φ28mm, Φ37mm
Application: Body sculpting, cellulite treatment, skin tightening, deep wrinkle removal, pore contraction, improved skin metabolism, and reduction of oily acne and under-eye puffiness.



Training and onboarding

Every system includes:

- Comprehensive training materials, including a detailed device operation and treatment guide
- Support with treatment protocols and recommended parameters
- Guidance on device handling, safety, and maintenance
- Ongoing technical support from the Power to Practice team

Warranty and maintenance

- 12-month manufacturer warranty included with all systems
- Servicing is available through Power to Practice as an additional support option
- Recommended servicing intervals vary based on usage, typically every 12 months

Why choose the HS-900?

This single system replaces the need for multiple devices — giving your clinic:

- A broader service offering
- Lower capital investment
- More space and workflow efficiency
- Consistent training and support from one supplier

It's not just a laser — it's a growth platform for your practice.

power to practice

Contact

info@powertopractice.com.au

www.powertopractice.com.au

